



# Dispersion Technology Inc.

Characterization of Concentrated Dispersions and Emulsions, Liquids and Porous materials

## Model DT-1202 Acoustic and Electroacoustic spectrometer: Particle sizing and Zeta potential measurement in concentrates.



### Wide Range of Applications:

- Nanotechnology
- Colloid stability
- Ceramic slurries
- Cement slurries
- Battery slurries
- CMP slurries
- Cosmetics
- Paints and pigments
- Non-aqueous systems
- Clays and minerals
- Food emulsions
- Mixed dispersions
- Structured systems
- Photo materials

**Model DT-1200** has two unique sensors: Acoustic and Electroacoustic.

**Acoustic sensor** characterizes *particle size distribution* by measuring ultrasound attenuation at set of frequencies from 1 to 100 MHz and sound speed. The same ultrasonic raw data can be used for calculating compressibility, elastic modulus and longitudinal viscosity of any liquid sample (when “Rheological option” is installed).

**Electroacoustic sensor** is built as a probe for measuring *ζ-potential* in concentrates without dilution. The same probe can be used for monitoring sedimentation kinetics and for characterizing porous materials (when “Porous materials option” is installed).

**These sensors** can function either separately for individual measurements, or together, providing certain synergism in sample characterization.

## Available Options:

**Titration option** with one or two burettes allows conducting of complicated experiments involving modification of chemical composition of the liquid medium according to a certain protocol. There are several different types of protocols available: “pH ramp”, “pH stat”, “surfactant titration”, “temperature titration”. Titration “pH ramp” allows scanning of a certain pH range in single or multiple sweeps and usually performed for determining iso-electric point. Titration “pH stat” monitors amount of a particular reagent that is required for maintaining given pH. Surfactant titration reflects changes in  $\zeta$ -potential, particle size distribution, or both, with incrementally increasing surfactant concentration. It is used for determining optimum surfactant dose. Temperature titration requires installation of “heating control option”, which would allow performing T sweeps within a range from room T up to 50 C.

**Conductivity aqueous option** allows for measuring electric conductivity of aqueous systems within a range from  $10^{-3}$  to 10 S/m. This probe functions at MHz range and, consequently, is not affected by electrode polarization. The same probe is used for measuring porosity of a porous material if Porous materials option is installed.

**Conductivity non-aqueous option** allows for measuring conductivity of various solvents including non-polar liquids within the range from  $10^{-11}$  up to  $10^{-4}$  S/m. This option is identical in function to the DT-700 model. This option requires installation of “non-aqueous media option”, which is important for protecting instrument sensor from aggressive solvents if they are intended to be used.

**Rheological option** allows calculation of high frequency (MHz range) longitudinal rheological parameters such as compressibility, elastic modulus, viscosity, and performs test on Newtonian nature of the liquid sample.

**Porous materials option** allows characterization of porosity using the aqueous conductivity probe, as well as pore size and zeta potential of a porous material with electroacoustic probe. Characterization of these last two parameters would require calibration.

**External pump option** is required when viscous samples are monitored continuously, which can serve as a laboratory prototype for on-line characterization.

### N o m i n a l   S p e c i f i c a t i o n s :

Calculated parameters		Sample volume, minimum [ml]	
Mean size [microns]	0.005-1000	Size only	15
Zeta potential [mV]	$\pm(0.5\%+0.1)$	Zeta only	0.1
Weight fraction / porosity	$\pm 0.1\%$	Both, no sedimentation	15 +0.1
Compressibility E10 [1/Pa]	$\pm 0.003$	Both with mixing	70
Bulk viscosity [cP]	$\pm 0.01$	Both with titration	110
Debye length [nm]	$\pm 0.1$	Both with pumping	150
Measured parameters		Sample requirements	
Temperature [C 0]	0 to 100, $\pm 0.1$	Volume fraction, % (1)	0.1-50
pH	0.5-13.5, $\pm 0.1$	Conductivity	none
Frequency range [MHz]	1-100	pH	0.5-13
Ultrasound attenuation [dB/cm/MHz]	0 to 20, $\pm 0.01$	Temperature [C 0]	<50
Sound speed [m/sec]	500 to 3000, $\pm 0.1$	Viscosity of media [cP] (2)	<20,000
Electroacoustic signal [mV(s/g) <sup>1/2</sup> ]	$\pm 1\%$	Viscosity of sample [cP]	<20,000
Conductivity [S/m]	$10^{-11}$ to 1, $\pm 1\%$	Particle size [microns]	0.005-1000
Measurement time [min]	0.5-10	Zeta potential [mV]	none

(1) Instrument can measure attenuation well above 50% vl, but verification of the theory for computing particle size and zeta potential is not possible above this limit.

(2) The “micro-viscosity” is important for theoretical calculation. It might be different than “macroscopic” viscosity for gels and other structured systems measured with conventional rheometers.

**Physical Specifications.** Electronic unit: weight 20 kG, sensor unit 30 kG. Power: 100-250 VAC, 50-60 Hz, <300 W.  
Software: embedded Windows HP, MS Office optional

Mail: 364 Adams Street, Bedford Hills, NY 10507  
Tel. 914-241-4777, Fax 914-241-4842

Email: [info@dispersion.com](mailto:info@dispersion.com)  
Homepage: [dispersion.com](http://dispersion.com)