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A figura 5 shows the particle size distributions of an active ingredient, domperidone, and excipients, aluminum and silicon oxides (alumina and silica) and magnesium carbonate


Figure 4: Particle size distributions for an active ingredient and excipients


Figure 6: Particle size distribution of a micro-emulsion

In the case shown in figure 5.a, a concentrated emulsion has been placed on a microscope plate. The specimen was carefully prepared in view of obtaining fully separated drops and ensuring optimal conditions for the analysis.

For this type of image, digital processing (figure 5.b) involves thresholding to convert the image into black and white. A selection filter was then applied to select isolated particles and exclude agglomerates. In the example shown, the filter applied allowed the most spherical particles to be selected.


Figure 5: Processing of an image of a microemulsion used in cosmetics

a) original image b) processed image

The size of the drops of this emulsion was determined on the basis of the analysis of 1049 particles. The corresponding particle size distribution is shown in figure 6.

The typical diameters of this monomode and monodisperse distribution are as follows:

$$D_{10} = 7.1 \mu\text{m}, D_{50} = 11.7 \mu\text{m}, D_{90} = 16.6 \mu\text{m}$$

## Conclusion

Whether handling emulsions, dry powders or powders dispersed in a liquid, diffraction laser particle size analysis combined with morphological analysis is an indispensable tool at all phases in the life of a drug.

Research and development phases, the request for marketing authorization or the acceptance inspection of raw materials or the inspection of the quality of the production of pharmaceutical preparations can be analyzed using these methods.

## References

- [1] International Standard , ISO 13320-1, "Particle Size Analysis - Laser Diffraction Methods", General Principles
- [2] International Standard , ISO 13322-1, "Particle Size Analysis - Image Analysis Methods Part 1 : Static Image Analysis Methods"
- [3] United State of Pharmacopeia : USP <429> Light Diffraction Measurement of Particle Size
- [4] United State of Pharmacopeia : USP <776> Optical Microscopy
- [5] United State of Pharmacopeia : USP <1058> Analytical Instrument Qualification
- [6] Intrinsic viscosity and the polarizability of particles having a wide range of shapes  
Advances in Chemical Physics, Volume XCI, (1995)
- [7] Compaction Simulator Studies of a New Drug Substance: Effect of Particle Size and Shape, and Its Binary Mixtures with Microcrystalline Cellulose  
Pharmaceutical Development and Technology, Volume 1, Issue 2, p. 119-126, (1996)
- [8] What is the role of particle morphology in pharmaceutical powder aerosols  
Expert Opinion on Drug Delivery, Volume 5, Number 8, pp. 909-914, (2008)
- [9] Divided solids characterization as part of pharmaceutical development  
STP Pharma pratiques, Volume 13, Number 4, p. 245-251, (2003)